## Réfs. : A032002 - A032003 - A032004

- NOTICE D'UTILISATION DU THERMOMÈTRE NUMÉRIQUE
- OPERATING INSTRUCTIONS FOR DIGITAL THERMOMETER
- BEDIENUNGSANLEITUNG FÜR DAS DIGITALTHERMOMETER
  - HANDLEIDING VOOR DE DIGITALE THERMOMETER
- INSTRUCCIONES DE FUNCIONAMIENTO DEL TERMÓMETRO DIGITAL
- MANUAL DE INSTRUÇÕES PARA O TERMÓMETRO DIGITAL
- ISTRUZIONI OPERATIVE PER IL TERMOMETRO DIGITALE
  - DIGITÁLNÍ ELEKTRONICKÝ TEPLOMĚŘ
  - A DIGITÁLIS HŐMÉRŐ HASZNÁLATI ÚTMUTATÓJA
  - NOTITA DE UTILIZARE A TERMOMETRULUI DIGITAL
  - INSTRUKCJA OBSŁUGI TERMOMETRU CYFROWEGO

Babymoov  
Parc Industriel des Gravanches  
16, rue Jacqueline Auriol  
63051 Clermont-Ferrand cedex 2  
France  
www.babymoov.com

**Hotline (only for France) :**  
**04 73 28 37 37**

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## EN IMPORTANT. READ INSTRUCTIONS BEFORE USE AND KEEP FOR FUTURE REFERENCE

**FEATURES**

- Temperature Range: 32 °C - 42.9 °C (90 °F - 109.9 °F)
- Accuracy: ± 0.1 °C (35 °C - 42 °C); ± 0.2 °C (< 35 °C)
- Digital LCD display
- Battery: 1.5 VDC (LR41) alkaline
- Power consumption: 0.15 mW
- Dimensions: 118 mm (L) x 25 mm (W) x 10 mm (H)

**HOW TO TAKE TEMPERATURE?**

- Press power switch, the screen shows the full segment 188.8 °E with a sound of "B", which indicates sufficient power and normal self-test routine as well as for measurement.
- Memory function: after the full segment, the last measured value with the letter "M" are displayed, then to show 37 °C / 98.6 °F.
- Release the button, then displays « Lo °C » or « Lo °F » and ready to measure. During measurement the symbol « °C » or « °F » will flash while the temperature is climbing. This will also happen when temperature increase remains below 0.1 °C / 0.2 °F for about 16 seconds. If temperature remains below 32 °C / 90 °F, the symbol « °C » or « °F » will flash for about 60 seconds. Then it will auto-off for saving power.
- A buzzing tone will be heard for about 10 seconds which means the finish of measurement. Normal alarm: sound of « bi...bi...bi...bi » means normal temperature that below 37,8 °C / 100 °F. Fever alarm: sound of « bi...bi...bi » means fever that the temperature is above 37,8 °C / 100 °F.
- During measurement LCD will display "Lo °C" or "Lo °F" if the measured temperature is below 32,0 °C (90 °F), or will display "Hi °C" or "Hi °F" if the measured temperature is above 43 °C (109,9 °F).
- Turn off the unit by just slightly pressing power switch, or leave it for 9 minutes, it will auto-off for saving power.
- Children should be under close supervision at all times while using this product. After monitoring the temperature, the unit should be removed from the child immediately.
- "°C" or "°F" exchange function: if the power switch was pushed more than 2 sec after power-on, the unit "°C" (°F) will be exchanged.

**DIFFERENT METHODS OF USE**

- You are advised to measure your temperature at specific time in morning, if it is earlier or later, the temperature would be higher or higher.
- **Used in mouth.** Please close your mouth for 1 minute for stability of temperature, breath calmly through nose. Respiratory air would affect the accuracy of measure of your mouth. Normal temperature: about 37 °C / 98,6 °F.
- **Used in anus of rectum.** This is the way that can get the most steady temperature and suits for both infants and children. The sensor shall be inserted for 2-3 cm. Vaseline can be used for easier inserting. Normal temperature: about 37,6 °C / 99,7 °F.
- **Used under armpit.** Clean the thermometer with a dry cloth and put it under armpit. Press and fix it with the arm. The temperature value of this measurement is not so accurate that no recommended for accurate measurement. Normal temperature: about 36,7 °C / 98,1 °F.

**Notes:** the normal temperatures are different from the different human bodies and measuring time. Sterilization is required before every measurement.

**BATTERY REPLACEMENT**

Battery must be replaced if  symbol appears during measurement. Make sure a battery of the same type is on hand and proceed as follows: Remove battery compartment lid and change with a button battery (LR41). Ensure battery is correctly positioned with "+" facing up.

**WARNING:** Keep battery out of children's reach. swallowing battery could be fatal. battery should not be charged or placed into extreme heat as it may explode.

**CLEANING AND MAINTENANCE**

- The thermometer is waterproof and can be cleaned by water or sterilization lotion with soft cloth; Sterilization lotion below is recommended:
- Name:** Density - Medical alcohol: 95 %
- Notes:** other lotion, density or cleaning method may damage the thermometer.

**CAUTIONS**

- Store the unit in the protective case while not in use.
- Do not store or keep in direct sunlight, or in high temperatures.
- Keep away from strong electromagnetic source to avoid interference.
- Do not open the thermometer except to replace battery.
- Do not drop the unit, as this may ruin the unit.
- It shall be operated in an environment at the temperature of 10 °C / 50 °F and the humidity of 40-80 % while stored and transported in an environment at the temperature of -10 °C - 40 °C (14 °F - 104 °F) and humidity of 30 - 90 %
- Never put it in boiled water for sterilizing.
- Please take out the battery if not use the unit for longtime.

**REMARKS**

- It is suitable for measuring temperature of human body and necessary of families.
- It is designed and manufactured complying with relevant requirements in ISO 80601-2-56.
- It has passed EMC emission test, immunity test according to EN60601-1-2 and biocompatibility test according to MDD annex I.
- According to the degree of protection against electric shock, it is classified to BF, in accordance with EN60601-1.
- Electrical and electronic items and electronic items are collected separately. Do not discard unwanted electrical and electronic items with unsorted municipal waste, but ensure that they are collected separately.
- Please replace battery before claim for replacing a new product.
- The scrap of product and battery are to be disposed according to relevant national regulations

**Caution!** Be sure that both the battery and the thermometer are kept away from children.

**Complies with safety requirements**

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